

Every OneTouch® meter is designed to make it easy to test your blood glucose and help you manage your diabetes.

Step
1

Insert a OneTouch Verio® test strip



Step
2

Apply blood to either side of the strip



Virtually pain-free testing

And that's it.

Your results appear on screen. Our ColorSure® technology will instantly show if your results are **in or out of range**.



Connect to the **OneTouch Reveal® App** for even more insight. See reverse side.

Before using the OneTouch® meter to test your blood glucose, carefully read the Owner's Booklet and inserts that come with the components of the system. Treatment decisions should be based on numerical result and healthcare professional's recommendation and not solely on where the result falls within the range limits.

Sync data, seamlessly, from your OneTouch® meter directly to your smartphone. Let's get started....

1. Download the FREE OneTouch Reveal® app

- Visit the app store on your mobile device, search "OneTouch Reveal", then install the app



2. Activate Bluetooth® on your OneTouch® meter

For the OneTouch Verio Reflect® meter

- Turn your meter on. Press and hold **OK** until the startup screen appears
- Go to the settings menu > select and activate Bluetooth®

For the OneTouch Verio Flex® meter

- Press the **OK** button to turn your meter on
- Press the **OK** and **^** buttons at the same time to activate Bluetooth®



Bluetooth®

Don't forget to also activate Bluetooth® on your mobile device.

3. Connect

- Open the OneTouch Reveal® app on your mobile device and follow the on-screen instructions
- When prompted, enter the PIN shown on your meter screen into your mobile device and tap **Pair**



You can also view your account on your computer at OneTouchReveal.com

12 3456

Example of a PIN number

Now you're ready to test and sync!

Discover more at OneTouch.com or call **1-866-693-0599**.