OneTouch Reveal[™] web app Report Reference Guide

Your guide to the different reports offered by the OneTouch Reveal™ web app for healthcare professionals







Verio Reflect® meter



Verio Flex® meter



Here to Help

In this guide, you'll find examples of some of the key OneTouch Reveal™ web app reports with detailed descriptions of the different information presented in each.

The reports display data in charts and graphs, highlight patterns and provide key statistics that can help reveal meaningful insights.

By showing them the bigger picture, the reports can encourage your patients to stay on track and help them to better manage their blood glucose between visits.



The OneTouch Reveal[™] web app helps to see more information that can drive therapy decisions.

Contents


Patient Summary	
Progress Report and 14-Day Logbook	2
Excursion Analysis	5
Logbook	6
Data List	7
Other Reports	8
Glossary	Ç



Patient Summary Report

Gives you an at-a-glance view of your patient's glycemic control over a selected date range.

- (1) Review prioritized pattern messages so you know what to focus on first.
- 2) Compare glycemic statistics from two separate time periods to see patient progress.
- 3 View percentages of pre- and post-meal glucose readings in and out of your patient's target ranges.
- 4) Review a time-of-day chart with a personalized patient schedule along with glucose monitoring and insulin dosing statistics.






Progress Report and 14-day Logbook

Optimize patient education by seeing what they see.

Shows you a version of the same report your patients can generate from their OneTouch Reveal™ mobile app.

- 1 Get at-a-glance view of your patient's glycemic control over a 14-day period and compare with previous period to see patient progress.
- 2 View patients blood glucose readings in a trendline.
- 3 See glucose readings along with tags and notes applied by the patient in the logbook.






Excursion Analysis Report

Summarizes patterns or variability above or below patient's target range.

- 1) See times of day when your patient is in and out of range.
- 2 Review excursions by date to track your patient's progress day by day.
- 3 See weekly patterns of low, high or variable blood glucose readings.







Logbook Report

Displays blood glucose readings, bolus doses and other events with columns as time of the day and rows as day of the week (or date) in familiar logbook format.

- See your patient's before and after meal blood glucose ranges.
- 2 See glucose readings organized by patient's personalized schedule and sorted chronologically by date.
- 3 Trace blood glucose reading from preceding events such as insulin bolus or carb intake.
- 4 Statistics include
 - Average glucose value
 - Total insulin dose
 - Total carbohydrate intake







Data List Report

Displays a tabular view of all data captured by the devices for a selected date range.

1) View an extensive list of data captured by connected devices e.g. high and low glucose results, insulin doses, etc. by time of day including meal tags applied by patient.

Sort by date order or result type.

ONETOUCH Reveal®

Data List: 11/11/2020 - 11/24/2020

Patient 4. Jack Patient ID: TTCP4 Overall Target: 70 - 180 mg/dL DOB: 3/1/1982 Clinic: Trident Training Clinic Before Meal Target: 70 - 130 mg/dL After Meal Target: Clinician: Clinic Manager, TTC Gender: Male 70 - 180 mg/dL

Blood sugar readings that you manually entered on the mobile app are not included in reports generated on this web app.

Warning: Reports are intended for use by healthcare professionals with experience treating patients with diabetes only.

Date	Time	Timeslot	Result Type	Value	Serial #	Comments
11/24/2020	5:22 pm	Before Dinner	Glucose	193 mg/dL	ZKP46FDBY	Before Meal
11/24/2020	1:40 pm	Before Lunch	Glucose	340 mg/dL	ZKP46FDBY	Before Meal
11/24/2020	11:44 am	Before Lunch	Glucose	188 mg/dL	ZKP46FDBY	Before Meal
11/24/2020	1:22 am	Overnight	Glucose	176 mg/dL	ZKP46FDBY	
11/23/2020	9:45 pm	After Dinner	Glucose	155 mg/dL	ZKP46FDBY	After Meal
11/23/2020	2:25 pm	After Lunch	Glucose	138 mg/dL	ZKP46FDBY	After Meal
11/23/2020	9:58 am	After Breakfast	Glucose	57 mg/dL	ZKP46FDBY	After Meal
11/23/2020	8:06 am	Before Breakfast	Glucose	110 mg/dL	ZKP46FDBY	Before Meal
11/23/2020	12:42 am	Overnight	Glucose	199 mg/dL	ZKP46FDBY	
11/23/2020	12:17 am	Overnight	Glucose	99 mg/dL	ZKP46FDBY	
11/22/2020	6:41 pm	Before Dinner	Glucose	76 mg/dL	ZKP46FDBY	Before Meal
11/22/2020	6:12 pm	Before Dinner	Glucose	129 mg/dL	ZKP46FDBY	Before Meal
11/22/2020	4:04 pm	After Lunch	Glucose	53 mg/dL	ZKP46FDBY	After Meal
11/22/2020	3:17 pm	After Lunch	Glucose	101 mg/dL	ZKP46FDBY	After Meal
11/22/2020	12:38 pm	Before Lunch	Glucose	262 mg/dL	ZKP46FDBY	Before Meal
11/22/2020	11:59 am	Before Lunch	Glucose	222 mg/dL	ZKP46FDBY	Before Meal
11/22/2020	9:22 am	After Breakfast	Glucose	96 mg/dL	ZKP46FDBY	After Meal
11/22/2020	8:39 am	Before Breakfast	Glucose	212 mg/dL	ZKP46FDBY	Before Meal
11/21/2020	11:01 pm	Bedtime	Glucose	158 mg/dL	ZKP46FDBY	
11/21/2020	11:00 pm	Bedtime	Glucose	273 mg/dL	ZKP46FDBY	
11/21/2020	8:20 pm	After Dinner	Glucose	70 mg/dL	ZKP46FDBY	After Meal
11/21/2020	2:25 pm	After Lunch	Glucose	147 mg/dL	ZKP46FDBY	After Meal
11/21/2020	10:16 am	After Breakfast	Glucose	148 mg/dL	ZKP46FDBY	After Meal
11/21/2020	1:54 am	Overnight	Glucose	49 mg/dL	ZKP46FDBY	
11/20/2020	8:31 pm	After Dinner	Glucose	172 mg/dL	ZKP46FDBY	After Meal
11/20/2020	7:16 pm	After Dinner	Glucose	88 mg/dL	ZKP46FDBY	After Meal
11/20/2020	2:06 pm	After Lunch	Glucose	76 mg/dL	ZKP46FDBY	After Meal
11/20/2020	12:08 pm	Before Lunch	Glucose	208 mg/dL	ZKP46FDBY	Before Meal
11/20/2020	6:43 am	Before Breakfast	Glucose	78 mg/dL	ZKP46FDBY	Before Meal
11/19/2020	9:11 pm	After Dinner	Glucose	369 mg/dL	ZKP46FDBY	After Meal
11/19/2020	2:20 pm	After Lunch	Glucose	85 mg/dL	ZKP46FDBY	After Meal
11/19/2020	9:49 am	After Breakfast	Glucose	109 mg/dL	ZKP46FDBY	After Meal
11/19/2020	5:48 am	Overnight	Glucose	65 mg/dL	ZKP46FDBY	
11/18/2020	7:31 pm	After Dinner	Glucose	121 mg/dL	ZKP46FDBY	After Meal
11/18/2020	11:48 am	Before Lunch	Glucose	217 mg/dL	ZKP46FDBY	Before Meal
11/18/2020	9:24 am	After Breakfast	Glucose	78 mg/dL	ZKP46FDBY	After Meal
11/18/2020	1:21 am	Overnight	Glucose	197 mg/dL	ZKP46FDBY	
11/17/2020	9:08 pm	After Dinner	Glucose	59 mg/dL	ZKP46FDBY	After Meal

Reported for Patient 4, Jack on 11/24/2020 from OneTouch Reveal® 4.3.1 OneTouch® Ultra®2 ZKP46FDBY

Page 1 of 2




Other Reports

The OneTouch Reveal™ web app also offers a number of other reports that you can use.

Integrated Daily View Report

Provides a time-of-day view of blood glucose readings, insulin and carbs



Details by Day of Week Report

Provides a summary of blood glucose, insulin and carbs by day of week



Details by Time of Day Report

Provides a summary of blood glucose, insulin and carbs by time of day





Glossary

Definitions of terms used in OneTouch Reveal[™] web app reports

Average daily risk range (ADRR). A measurement of variability in blood glucose fluctuations, using a range scale of 0-60. A measurement of 0-19 represents low risk, 20-39 represents medium risk, and 40-60represents high risk.

Basal dose. The continuous amount of insulin the body needs throughout the day.

BG. Abbreviation for blood glucose.

Bolus dose. The additional units of insulin needed to cover carbohydrates or to correct a high blood glucose level.

Coefficient of variation (CV). Defined as the ratio of the standard deviation (SD) to the mean (average), expressed as a percentage: %CV = 100 X SD/mean.

Data List Report. A report that provides a tabular view of all data captured by the device for a selected date range, displayed in chronological order.

Data record. A unit of information in OneTouch Reveal[™] web app such as a blood glucose reading, an insulin dose, or other information (e.g. a name).

Date range. The number of days of data (14, 30, 90, or a custom range) contained in each report.

Details by Day of Week Report. A report that provides a summary of glucose readings, insulin, and carbohydrate intake by day of week, for a selected date range.

Details by Time of Day Report. A report that provides a summary of glucose readings, insulin, and carbohydrate intake by time of day, for a selected date range.

Device. Blood glucose meter or insulin pump that can transfer data to the OneTouch Reveal™ web app

Excursion Analysis Report. A report that provides a listing of hyperglycemia, hypoglycemia, variability, and pump patterns for a selected day range.

Footnotes. Messages appearing at the bottom of reports that communicate additional information.

Glucose average. The arithmetic mean calculated for a set of glucose readings.

High blood glucose indicator (HBGI). A measure for estimating hyperglycemia risk. A score of 4.5 or lower indicates low risk, a score between 4.5 and 9.0 indicates moderate risk, and a score higher than 9.0 indicates high risk.

HIGH blood glucose readings. When 'HIGH' (or 'HI') appears in a report in capital letters, it refers to glucose readings higher than the measurement range of the meter. Consult the Owner's Booklet for the measurement range of your device.

Integrated Daily View Report. A report that provides glucose readings daily CGM tracings with basal, I:C (insulin to carbohydrate) ratio, ISF (insulin sensitivity factor), total daily carbohydrates, and insulin dose data by time of day, for a selected day range.

Logbook Report. A report that lists glucose readings and pump, carbohydrate, and insulin data in columns by time of day, for a selected date range.

Low blood glucose indicator (LBGI). A measure for estimating hypoglycemia risk. A score of 1.1 or lower indicates minimal risk, a score between 1.1 and 2.5 indicates low risk, a score between 2.5 and 5.0 indicates moderate risk, and a score higher than 5.0 indicates high risk.

Low blood glucose readings. When 'LOW' (or 'LO') appears in a report in capital letters, it refers to glucose readings lower than the measurement range of the meter. Consult the Owner's Booklet for the measurement range of vour device.

Median. The middle value in a data set (taken as the average of the two middle values when the sequence has an even number of values).

Meter ID. A serial number stored in the memory of each meter. OneTouch Reveal[™] web app retrieves the meter ID when transferring meter results so it can keep track of the meter from which the data originated.

Patient Summary Report. A report that provides an overview of glucose patterns and excursions, testing and dosing regimens, and key comparative statistics for a selected date range.

Pattern messages. Messages appearing in reports that may help identify trends in patient data.

Schedule. A 24-hour day divided into eight time periods that can be customised to meet an individual's personal daily routine.

Standard deviation (SD). A measure of dispersion – i.e. how much the test results in a certain set are scattered around the mean. A low SD signifies that the test results are tightly clustered; a high SD signifies the results are widely scattered.

Tags. A note attached to a result to further identify the data. The tag may indicate that the result is a Before Meal, After Meal, Fasting, or Bedtime result.

Target range. The range (upper and lower limits) of preferred glucose levels.

Time period. The eight periods within a 24-hour day used to organize data transferred to the OneTouch Reveal[™] web app.

Transfer. The procedure that moves data from a meter or insulin pump to the OneTouch Reveal[™] web app.

Unit of measure. Blood glucose readings as well as other test results are reported in mg/dL or mmol/L.



Discover the benefits of the OneTouch Reveal™ app for you and your patients



OneTouch Verio Reflect® and OneTouch Verio Flex® meters

connect to a single tool that downloads and seamlessly aggregates data.

OneTouch Reveal[™] mobile app.

Changes the way your patients see their blood sugar.

- Effortlessly transforms results into quick color snapshots
- Identifies high, low and recurring patterns and
- Automatically sends notifications

OneTouch Reveal[™] web app. A simple way to help patients stay on track — in the moment and between visits:

- Simple, colorful visuals that highlight patterns
- Connects blood glucose, food and insulin in new ways
- Results can be shared with you, between visits or during appointments



